

VarioVetPro - electrically height adjustable X-ray, treatment and surgery table

→ Especially developed for practices with limited space

Brief description

- Serves as X-ray, treatment and surgery table
- Electrically height adjustable table with two columns
- Working height: 80 to 120 cm
- Integrated tilting mechanism allows to lift the table top on both sides via two switches
- The cassette / grid drawer under the table top (for unit sizes up to 35 x 43 cm) provides for the use as X-ray table
- Table top is made of white, synthetic material, which is easy to clean and disinfectant proof
- Dimensions of table top: 144 x 65 cm
- Table is moveable on wheels and allows an easy relocation within the veterinary practice



Integrated tilting mechanism

With either hand- or foot switch



VarioVetPro - electrically height adjustable X-ray, treatment and surgery table



Cassette drawer for
cassette sizes up to
35 x 43 cm

Hand or foot switches
for electrical adjustment
of the table top /
working height

Specifications subject to revision without notice

Specifications subject to revision without notice
The editor strives to impart correct and up to date information. The provided specifications are based on current knowledge and are subject to revision without notice. This brochure is subject to correction. The editor assumes no responsibility for the information being up to date, correct and complete.

All furnished logos, pictures and graphics are property of the particular company and subject to copyright of the licensor. Use, dissemination, distribution or copying of the pictures, logos or text compiled or processed by the editor is subject to our written consent. All rights reserved.

Distributor